

# A Noninvasive Prenatal Screen that Achieves $\geq 4\%$ Fetal Fraction in $>99.9\%$ of Patients

Rotem Ben-Shachar, PhD; Katie Johansen Taber, PhD; Dale Muzzey, PhD; Susan Hancock, MS, CGC

All authors were employed by Myriad Genetics, Inc. at the time of this study

## PURPOSE

- For millions of pregnant women, noninvasive prenatal screening (NIPS) based on cell-free DNA (cfDNA) detects whether their pregnancies are at elevated risk for fetal chromosomal abnormalities.
- Fetal fraction (FF), the proportion of cfDNA originating from the placenta, can impact the accuracy of NIPS, and many laboratories fail samples with low FF, commonly defined as FF  $<4\%$ .
- FF has been shown to negatively correlate with body mass index (BMI), pregnancies with trisomy 18 or 13, and early gestational age, resulting in higher test failure rates in these populations.
- A whole-genome sequencing (WGS)-based NIPS that employs FF amplification (FFA) technology for all samples has been shown to increase FF by 3.9-fold for samples with low FF.


## METHODS

- We retrospectively analyzed results from 104,710 patients who underwent NIPS with FFA during an eight-month period.
- The FFA technology increased FF by preferentially sequencing short cfDNA fragments, known to be enriched for fetal-derived cfDNA.
- FF was assessed for patients who received a screening result (N= 104,557).
- BMI data were available for 65,773 patients.

## RESULTS

Figure 1. FFA performance in  $>100,000$  clinical samples.

(A) Fetal fraction distribution across  $>100,000$  clinical samples. (B) FF was assessed looking at the 1st percentile FF for a variety of combinations of BMI and gestational age. The dot size corresponds to the number of women with the indicated BMI and gestational age.



- Median maternal age was 31 years and median gestational age was 12 weeks.
- No patients had tests failed due to FF  $<4\%$ .
- Less than 0.1% of patients had FF below 4%.
- Ninety-nine percent of patients had FF  $>6.4\%$  (Fig 1A, dashed line).
- In patients with gestational age less than 12 weeks (N= 48,722), 99% of patients had FF  $>6\%$ .
- In patients with BMI  $\geq 30$  (N=23,380), 99% of patients had FF  $>5.4\%$ .
- Even in patients with BMI  $\geq 40$  and gestational age of 10 weeks, 99% of patients had FF  $>4.6\%$ .

## CONCLUSIONS

- A commercial NIPS using high-throughput FFA achieves sufficiently high FF levels to provide confident results regardless of a woman's risk factors, preventing unnecessary test failure.