

# Early Outcomes from pcfDNA Screening at 8-10 weeks Gestation for Samples with Additional Findings

Sarah Hash<sup>1</sup> and Carly Hammer<sup>1</sup>

<sup>1</sup> Myriad Genetics, Inc., Salt Lake City, UT



## Background

- Historically, prenatal cell-free DNA (pcfDNA) screening has been available starting at 10 weeks gestation
- Recent advancements have enabled screening as early as 8 weeks gestation (8w0d–9w6d)
- This study evaluates early pregnancy outcomes associated with “Additional Findings” results in pcfDNA samples collected during this earlier gestational window

## Methods

- Data from January 2025 to October 2025 were analyzed to identify trends in early pregnancy outcomes
- Laboratory genetic counselors contacted ordering providers for all pcfDNA samples collected between 8w0d–9w6d gestation that resulted with “Additional Findings”:
  - Suspected interfering maternal findings
  - Rare autosomal trisomies (RATs; trisomies in autosomes other than 21, 18, and 13)
  - Other complex findings
- Genetic counselors documented:
  - Panel utilization
  - Details of the abnormal result and gestational age at collection
  - Pregnancy viability at the time of sample collection
  - Viability at follow-up phone call (~10–11 weeks gestation)
  - Presence of a reported vanishing twin


**Table 1.** Additional Findings in case set

Maternal Chromosome Findings	17
Maternal Mosaicism for XO	8
Maternal Mosaicism for XXX	2
Maternal chromosome 22 duplication	6
Maternal chromosome 22 deletion	1
Multiple Chromosome Findings	18
Multiple Trisomies	5
Multiple Monosomies	0
Monosomies and Trisomies	6
Partial deletion and duplication suggestive of a translocation	7
Singular Chromosome Finding	85
Single Trisomy	73
Single Monosomy	0
Mosaic or Ring X chromosome	3
Single Deletion	5
Single Duplication	4


## Results

- A total of 120 early gestational age (EGA) pregnancies with “Additional Findings” were identified (**Table 1**).
  - 17 cases (14%) resulted with maternal chromosomal findings, most commonly maternal mosaicism for monosomy X or a maternal chromosome 22q duplication
  - 18 cases (15%) showed multiple chromosomal abnormalities, most frequently involving partial deletions and duplications suggestive of a derivative chromosome inherited from a parent with a reciprocal translocation
  - 85 cases (71%) revealed a single chromosomal abnormality (**Figure 1**).
    - 73 identified as RATs, with the most prevalent being:
      - Trisomy 15 (21/73)
      - Trisomy 22 (13/73)
  - Optional panel utilization varied. Expanded Aneuploidy Analysis (EAA) was included in 68% of orders (82/120). Nearly half of individuals (48%, 57/120) received all three optional panels; Sex Chromosome Analysis (SCA), EAA, and Microdeletions (**Figure 2**).

**Figure 1.** Most prevalent whole chromosome RATs identified




**Figure 2.** Distribution of optional add-on panel combinations (N = 120)




All ordered panels automatically include the core chromosomes 21, 18, and 13.

- Excluding those results with maternal findings, 103 cases with Additional Findings were assessed for early pregnancy outcomes (**Figure 3**).
  - 66% (68/103) pregnancies viable at collection
  - 12% (12/103) pregnancies not viable at collection
  - 22% (23/103) with unknown viability at collection
- Follow-up data showed that of the 68 initially viable pregnancies:
  - 54% (37/68) remained viable
  - 41% (28/68) resulted in miscarriage
  - 4% (3/68) viable with a vanishing twin identified via ultrasound
- Among the 80 pregnancies with known follow-up data, 40 (50%) experienced miscarriage between 8- and 11-weeks' gestation (**Figure 4**).


**Figure 3.** Early outcomes from pcfDNA screening at 8-10 weeks gestation for samples with Additional Findings



**Figure 4.** Early pregnancy outcome data for EGA samples with Additional Findings



**Figure 5.** Rate of pregnancy loss at outcome collection by chromosome finding



## Conclusions

- pcfDNA screening at EGA provides clinically meaningful information for patients and providers
- Analysis of all chromosomal aneuploidies, rather than only the common trisomies, in pcfDNA screening performed at EGA may provide insight into the causes of early pregnancy loss
- Abnormal findings may:
  - Offer insight into potential etiologies of pregnancy loss
  - Facilitate targeted genetic counseling
  - Guide follow-up testing for both the pregnancy and the parents
  - Inform recurrence risk assessment for future pregnancies
- Many pregnancies with abnormal results experienced miscarriage before the gestational age at which pcfDNA screening has traditionally been available, identifying a population that previously could not benefit from this technology
- Utilization of pcfDNA screening between 8w0d-9w6d gestation can complement traditional diagnostic pathways to improve understanding of the genetic landscape of early pregnancy loss

**Disclosures:** All authors were employees of Myriad Genetics, Inc. at the time of this study and received salary and stock as compensation.

Email for questions: [Sarah.Hash@myriad.com](mailto:Sarah.Hash@myriad.com)