

# RNA Analysis from Residual Blood Aids the Interpretation of VUS Identified in Individuals Undergoing Hereditary Cancer Genetic testing




Paola Nix<sup>1</sup>; Erin Mundt<sup>1</sup>; Benjamin Roa<sup>1</sup>

<sup>1</sup> Myriad Genetics, Inc., Salt Lake City, UT

## Background

- Functional RNA studies are a useful tool to help resolve the classification of variants of uncertain significance (VUS) suspected to have an impact on mRNA splicing.



- Splicing mutations are a possible cause of gene disruption
- Functional RNA studies help determine the level of aberrant splicing that a variant might cause

- RNA studies may also help establish the phase of co-occurring variants in cases where clinical management differs depending on whether variants occur *in cis* or *in trans*.

- We recently validated the use of blood samples collected in EDTA tubes for RNA studies.
- Residual blood submitted for hereditary cancer (HC) DNA testing can yield sufficient material for RNA isolation and functional analysis by RT-PCR and cDNA sequencing.

## Methods

- RNA analysis was performed on residual EDTA blood samples submitted for HC panel testing that carried a heterozygous VUS predicted to impact splicing by in silico models or carried co-occurring variants where phase could not be determined by NGS sequence or haplotype review.
- Total RNA was isolated and subjected to RT-PCR and cDNA sequencing.
- Allele-specific transcripts were reviewed to identify and quantify RNA isoforms produced by either the wild-type or the variant allele.
- These data were in turn compared to RNA from control blood and relevant tissue samples.

## Results


- We performed RNA analysis on 56 selected cases between May 2024-April 2025.
- RNA and cDNA of sufficient quantity and quality were isolated from 41/56 (73%) residual EDTA blood samples.
- 38 samples were analyzed for splicing effects and 3 samples were analyzed for variant phase for a total of 41 samples tested.

- Of these, 22/38 (58%) demonstrated a complete splice defect and were upgraded from VUS to likely pathogenic, while 16/38 (42%) showed incomplete or inconclusive splice defects and retained a VUS classification.


- Three samples each carried two variants in *MUTYH*, *MLH1*, and *CHEK2*, respectively.

- RNA analysis determined the cis or trans orientation in each case, ultimately providing more specific clinical management to the patient.



**Figure 1.** Allele-specific RNA studies from residual blood help interpret VUS identified in a patient's germline DNA test



**Example 1:** *ATM* c.3078-5T>C



**Figure 2.** Confirming variant phase provides more precise cancer risk estimate



## Conclusions

- Under certain conditions, residual samples can be used for RNA analysis, thereby avoiding the burden of secondary sample submission.
- This pathway facilitates both qualitative and quantitative analysis of allele-specific transcripts.
- The addition of residual blood RNA studies has resulted in a 5-fold increase in the number of variants analyzed in our laboratory during a 12-month period.
- These functional studies contribute substantially to the reclassification of splicing variants, providing definitive test results to more patients in a streamlined manner.