

F E B R U A R Y 1 2 , 2 0 2 6

pcfDNA Screening with Fetal Fraction Amplification Between 8-10 Weeks: Laboratory Experience with >28k Samples

Ronit Lebor¹, Claire Miller¹, Jack Roach¹, Jhett Bordwell¹, Seyedmehdi Shojaee¹,
Cathy Wicklund¹, Summer Pierson¹, Katie Johansen Taber¹, Dale Muzzey¹

1. Myriad Genetics, Inc.

Disclosures

All authors were employees of Myriad Genetics, Inc. at the time of this study and received salary and stock as compensation

Objective

- Explain benefits of fetal fraction amplification (FFA) in prenatal cell-free DNA (pcfDNA) screening
- Review early laboratory experience with pcfDNA screening with FFA between 8w0d and 9w6d gestation


AMPLIFY™ technology increases fetal fraction by leveraging differences in cfDNA length

Molecules sequenced in traditional cfDNA

(Avg. length ~ 165nt)



FF = 15%




Molecules sequenced in cfDNA with FFA



(Avg. length ~ 140nt)



FF = 35%



Background: Early Gestational Age (EGA) Validation



Methods

- pcfDNA screens ordered between 11/19/2024 and 07/21/2025 considered for analysis
- Number of orders received, screen-positive results, and screen failures calculated
- Fetal fraction (FF) medians and first and third quartiles (Q1, Q3) for the 8- and 9-week gestational age period calculated
- BMI information analyzed when provided

Results of weeks 8 and 9

- 28,446 orders were placed for patients between 8w0d and 9w6d

Result type	702 screen positive reports (2.47% screen positive rate)
Aneuploidies	
Trisomy 21	140 (19.94%)
Trisomy 13	48 (6.84%)
Trisomy 18	39 (5.55%)
Sex chromosome aneuploidy (SCA)	234 (33.33%)
Rare autosomal aneuploidy (RAA)	292 (41.60%)
Microdeletions	
1p36	8 (1.14%)
22q11.2	4 (0.57%)
15q11.2	3 (0.43%)
5p	4 (0.57%)
4p	2 (0.29%)

*25 were screen positive for two or more abnormalities


Screen failures

Reason for failure	61 total no call reports (0.2% screen failure rate)
Insufficient FF	19 Samples (0.06%)
Other QC metrics	42 Samples (0.14%)

Week	Rate
Week 8	37 Samples (0.13%)
Week 9	24 Samples (0.08%)


Median FFs

- **Across the cohort:**
15.4% (Q1=11.4%; Q3=20.2%)
- **8th week:**
13.8% (Q1=10.3%; Q3=18.2%)
- **9th week:**
17.1% (Q1=12.9%; Q3=22.0%)



Fetal Fraction distribution at high BMI

- BMI ≥ 30 : 35.5% (n = 6,246)
- Median FF 12.5%
(Q1=9.3%; Q3=16.1%)



Clinical Utility

- Screening at earlier gestational ages facilitates genetic counseling and pregnancy management, such as earlier diagnostic testing
- pcfDNA screening at 8 weeks aligns with initial prenatal visits that commonly occur at approximately 8 weeks gestation, avoiding the need for a later blood draw
- Earlier screening has the potential to reduce anxiety and uncertainty earlier in pregnancy

Conclusions

- Laboratory experience with pcfDNA screening for pregnancies between 8w0d-9w6d weeks gestation:
 - Demonstrated sufficient FF for analysis
 - Maintained a low screen-failure rate across a broad population, including patients with high BMI
- Screening for fetal aneuploidy at an earlier gestational age allows patients and providers to gain clinical insights earlier, leading to earlier decision making

Thank you

Questions?